

F334 Distribution of Genes Coding for Aminoglycoside
Acetyltransferases among Gentamicin Resistance Bacteria Isolated
from Aquatic Environments

Hyo Shim Han*, Young Jong Lee¹, and Jae Sung Jung
Department of Biology, Sunchon National University,
¹Department of Clinical Pathology, Kwangju Health College

Most resistance to aminoglycosides result from the production of plasmid-encoded aminoglycoside modifying enzymes namely, acetyltransferases, nucleotidyltransferases and phosphotransferases. Studies on the distribution of antibiotic resistance genes will supply data not only of clinical and epidemiological implications but also information about gene transfer processes in natural environments. DNA probes are useful tools studying epidemiology of antimicrobial resistance gene. Recently, PCR methods are used for the identification of aminoglycoside modifying enzymes. In this study we report the distribution of genes for aminoglycoside acetyltransferase among the gentamicin resistance bacteria isolated from aquatic environments.

F801 The *YDR1/BUR6* heterodimeric complex is a
Transcriptional Repressor in Yeast

* Hyun-Sook Yoo, Sung-Chul, Oh and Jong-Gil Na
Department of Biology, Kunsan National University

The human heterodimeric Dr1/DRAP1 is a general transcription repressor extensively studied *in vitro*. The two proteins repress transcription by binding TATA-binding protein (TBP), inhibiting the formation of transcription preinitiation complex. Recently the counterparts of the human Dr1 and DRAP1 were identified and named *YDR1* and *BUR6*, respectively, in the yeast *Saccharomyces cerevisiae*. The *YDR1* gene is essential for cell viability. *BUR6* is also essential for cell viability. The human Dr1 gene can replace the *YDR1* gene *in vivo*. *YDR1* overexpression *in vivo* confers slow growth phenotype. The *YDR1* and *BUR6* can be copurified in a gel filtration column, using affinity-purified polyclonal antibody against both yDr1 and Bur6 polypeptides. The gel mobility shift assay showed that Ydr1 and Bur6 repress *in vitro* transcription and this repression of basal transcription can be overcome by TBP, a global regulator. These results indicate that *YDR1/BUR6* functions as a global repressor of transcription *in vivo* and directly targets TBP. We are currently mutagenizing *YDR1* and *BUR6* to study structure-function of the two genes.