## Purification and refolding of human erythropoietin receptor produced in *Escherichia coli*

<u>Tae Wan Kim</u><sup>1</sup>, Taek Jin Kang<sup>1</sup>, Jin Ho Ahn<sup>2</sup>, Nam Young Kim<sup>1</sup>, In Seok Oh<sup>1</sup>, Yeo Joon Yoon<sup>3</sup>, Cha Yong Choi<sup>1,2</sup> \*

School of Chemical Engineering, College of Engineering, Seoul National University<sup>1</sup>;

Interdisciplinary Program for Biochemical Engineering and Biotechnology,

College of Engineering, Seoul National University<sup>2</sup>;

School of chemical engineering and bioengineering, University of Ulsan<sup>3</sup>

TEL: 02-880-7528, FAX: 02-874-1206

## Abstract

Erythropoieton receptor (EpoR) is a member of the cytokine receptor super-family and is expressed primarily on hematopoietic cells. EpoR plays an important role in regulating red blood cell production. Due to its low expression level in human tissue, studies on the biochemical and biophysical properties of EpoR have been limited. In the present study, an extracellular domain of recombinant human EpoR (rh-EpoR-ED) was expressed in *Escherichia coli* as inclusion bodies. Purification was achieved by immobilized metal affinity chromatography (IMAC) under denaturing condition and purified receptor was refolded by solid phase refolding method. We obtained sufficient quantities of rh-EpoR-ED, which can be available for biochemical and biophysical characterization through our production and purification method.

## References

- 1. Terry Lappin (2003), The cellular biology of erythropoietin receptor, *The oncologist* **8**(suppl 1), 15-18.
- 2. Kevin W. Harris (1992), Ligand binding properties of the human erythropoietin receptor extracellular domain expressed in *Escherichia coli*, *J. biological chemistr*. **267**, 15205-15209.
- 3. H Rogl (1998), Refolding of *Escherichia coli* produced membrane protein inclusion bodies immobilised by nickel chelating chromatography, *FEBS Letters* **432**, 21-26.
- 4. Middelberg A. P. (2002), Preparative protein refolding, *Trends in Biotechnology* **20**(10), 437-443.