

P95

Identification and Characterization of Imipenem-Resistant Clinical Isolates Producing Carbapenemase

Chul Park^{1,2}, Keun Sik Baik², Hyung Rak Kim^{1,2}, Seong Chan Park² and Chi Nam Seong^{2*}

¹*Department of laboratory medicine, Saint Garlo medical Center, Suncheon 540-718, Republic of Korea*

²*Department of Biology, Suncheon National University, Suncheon 540-742, Republic of Korea*

Imipenem-resistant bacteria were isolated from clinical specimens taken from hospitalized patients in Suncheon, Korea. Fifty-six isolates were phylogenetically analyzed based on 16S rRNA gene and *gyrB* gene sequence comparisons. Strains were affiliated with *Pseudomonas aeruginosa* (29 strains; 51.8%), *Acinetobacter baumannii* (22 39.3%), *Pseudomonas monteilii* (1), *Enterobacter hormaechei* (1), *Pseudomonas beteli* (1), *Enterobacter cloacae* (1) and *Pseudomonas geniculata* (1). Antibiotic resistance phenotype and MIC of the isolates was analyzed with the disc diffusion method and Vitek card. Resistance genotype was analyzed via PCR amplification of β -lactamase gene and repetitive sequence based PCR using ERIC primer. Six genotypes, in which two or more types in one strain, were found; genotype D (35 strains), E (29), B (19), A (17), C (11) and F (6).

P96

High-Throughput Screening for the Identification of Marine Microorganism Secreting Agarase and Cellulase

Duwoon Kim, Jin-Kyung Song, Sook-Hyun Ji, Na-Young Kang, Hea-Jin Woo, Bok-Mi Jung, Keun Sik Baik¹, Sung-Joo Jung, Tai-Sun Shin and Chi Nam Seong^{1*}

Division of Food Science and Aqualife Medicine, Chonnam National University

¹*Department of Biology, Suncheon National University*

High-throughput screening was performed for the identification of marine bacteria secreting agarase and cellulase. The highest agarase and cellulase activity was found in the culture supernatant of *Pseudoalteromonas* strain. The optimal pH and temperature of agarase and cellulase were determined by a reducing sugar assay, where the optimal pH of two enzymes was 8 and the optimal temperatures were 40 and 50°C, respectively. A thin layer chromatography (TLC) was used for the determination of agarose and cellulose hydrolysis reaction after incubating 200 μ l of substrates (0.1% agarose and 0.2% carboxymethyl cellulose) and 25 μ l of culture supernatants (2.29 mg/ml) at 40 and 50°C for 2hr. The presence of low molecular size of digested cellulose and agarose on the TLC plate indicated that *Pseudoalteromonas* strain secretes significant quantities of agarase and cellulase into extracellular environment. The assay has the advantage that agarase and cellulase activity can be directly measured using liquid cultures grown in a microtiter plate instead of separation or purification steps and is fast and easy to perform more adaptable for screening of a large number of samples.